



OC6000 Oxygen Concentrator with E-tank Back-Up

Works With Any Existing Anesthesia System!

Supera Anesthesia Innovations Anesthesia Innovations is proud to bring the revolutionary Pureline® OC6000 O2 concentrator to the biomedical research community for the first time. The revolutionary system provides on-demand oxygen to any existing anesthesia system. While similar machines have been in use in Europe for years, Supera Anesthesia Innovations has perfected the technology, added safety features, and now is making this money-saving equipment available here. The built in E-tank manifold allows mobile use of the machine without power and protects your subject in case of a power failure. The Pureline® OC6000 O2 Concentrator is FDA approved, UL tested, and backed by the industry's best three year warranty.

Product Information

- Ready for use with and Supera Anesthesia Innovations Anesthesia System or your existing anesthesia equipment
- Complete oxygen concentrator, frame and valve assembly
- Extracts Medical grade oxygen (USP93%) from the air and delivers on demand- up to 5 LPM
- Uses the industry standard DISS male thread connection
- Safety back-up oxygen supply (standard E tank) with fully integrated O2 flush with industry standard 30 LPM maximum flow
- Ability to use your oxygen flush with a simple built in switch
- FDA approved, UL tested
- Oxygen purity and flow safety alarms
- Quiet 40-dBA operation
- Made in the USA, state-of-the-art manufacturing with ISO quality and TUV SUV safety certifications
- Sturdy A-frame base with compact 18" square footprint

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com