

Metabolic Cage

THE INDUSTRY'S FIRST DISPOSABLE METABOLIC CAGE THAT
DRAMATICALLY SIMPLIFIES STORAGE AND ASSEMBLY.

MAIN FEATURES

- » Specifically designed for mice
- » Arrives irradiated and ready to use
- » Made from 100% PET BPA-Free plastic
- » Single-use components eliminate variability due to degradation of the collection surface from repeated washings
- » Single-use hydrophobic collection element assures complete sample collection
- » Unit can stand alone or be docked in an Innovive rack
- » Uses standard size 15ml conical tubes for collection*

FACILITY BENEFITS

- » Utilizes less storage space
- » Easy to assemble
- » Fewer pieces to clean
- » Recyclable, single-use components reduce contamination
- » Designed for convenience and reproducibility

COMPONENTS

- » Lid with a large Manniweb® filter media and metal grate
- » Innocage® Mouse cage bottom
- » 300ml Aquavive® water bottle
- » Sipper tube and bottle seal system
- » Drip receptacle for water
- » Crumb receptacle for food
- » Food & bottle wire bar hopper
- » Wire bar platform
- » Support cone with hydrophobic single-use liner
- » Collection separator
- » Stand

* 15ml conical tubes sold separately



Image 1. Side view of an assembled Metabolic cage on a stand.



Image 2. Front view of an assembled Metabolic cage.



Image 3. Back view of an assembled Metabolic cage.



V09/19